



Be Right™



Glass fiber filter, pore size 1.5 µm, 47 mm diameter, 100 pieces

Produktnr.: 253000
DKK Pris: Kontakt os
Tilgængelig

for TSS (Total Suspended Solids)

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 µm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

Specifikationer

| | |
|---------------------------|------------------|
| Diameter: | 47 mm |
| Materiale: | Glass Microfiber |
| porestørrelse: | 1.5 µm |
| Volumen / pakkestørrelse: | 100 /pk |