



Be Right™



Acid Reagent Powder Pillows (pk/100)

Produktnr.:

104299

Farligt

DKK Pris:

Kontakt os

Tilgængelig



Powder reagent in sealed packages for the determination of high range Silica, Silicomolybdate Method. Also used for the determination of high range Hypchlorite, Iodometric Method.

Specifikationer

Antal test:	100
Metodenavn:	Silicomolybdate
Metodenummer:	8185
Parameter:	Silica
platform:	Pulverpuder
Volumen / pakkestørrelse:	100 /pk